

# Model ER825 C<sup>4</sup>D Detector Multi-Channel Contactless Conductivity Detector



Capacitively-coupled contactless conductivity detection (C<sup>4</sup>D) systems apply an AC waveform to a transmitter electrode around a tube (or adjacent to a microfluidic channel) in which electrophoretic or chromatographic flow is occurring. The AC signal capacitively couples into the sample, which conducts the signal to a nearby receiver electrode. This, much attenuated, AC signal is demodulated to provide a DC signal which is dependent on the conductivity of the sample between the electrodes. This process occurs inside a headstage module, which outputs a signal to the C<sup>4</sup>D Detector. The Detector removes signal offset, amplifies and filters the signal, and also generates the excitation waveform that is sent to the headstage.

The CE or chromatography tubing is passed through the C^4D Headstage where the AC signal is applied and received.

# Specifications

Channels:	2 — 8 (specify when ordering)
Input connectors:	10 pin LEMO compatible socket
Signal gain:	x1, x10, x100
Analog signal resolution:	16 bits (Chart software) 24 bits (PowerChrom software)
Virtual serial port resolution	24 bit
Input gain ranges:	±20, 50, 100, 200, 500 mV ±1, 2, 5 V
Excitation frequency:	20 — 2000 kHz
Excitation amplitude:	2 — 40 V AC p-p, sinusoidal, @ <50 mA
Offset:	±5.5 V on ±10 V output range

WARRANTY: eDAQ Hardware units are supported by a three year warranty

#### www.**eDAQ**.com

E-mail: info@edaq.com

e-corder is a registered trademark of eDAQ Pty Ltd. All other trademarks are the property of their respective owners.

### Applications

- Capillary electrophoresis (CE)
- Ion chromatography (IC)
- Salt gradient chromatography
- Packing consistency for monolithic columns

## Compatibility

This Detector has a virtual serial interface (which that runs over its USB connection to the computer). Settings can be optimized with supplied C<sup>4</sup>D Profiler V2 software. The signal can be collected from the analog signal output, or via the USB/virtual serial interface. Many commercial data acquisition systems and chromatography recording systems will be compatible.

The ER815 is compatible with most CE and IC systems, including Agilent, Beckman Coulter and PrinCE instruments.

Compatible with all eDAQ headstage and platforms including:

- ET120 C4D Headstage for 360 365 µm OD capillary
- $\bullet$  ET125 C4D Headstage for chromatography tubing (1/16 inch OD or less)
- ET121 Microfluidic Platform for microchip electrophoresis.

Low pass filters:	1000 (off),10, 1 Hz
Output signals:	$\pm 5$ V max, or USB with virtual serial interface
Back panel connectors:	BNC (analog voltage) USB (virtual RS232 serial)
Power requirements: (mains adaptor supplied)	12 V DC, ~10 W
Dimensions (h × w × d):	65 × 200 mm × 250mm (2.6 × 7.9 × 9.8")
Weight:	1.8 kg (4.3 lb) maximum
Operating temperature:	0 to 35 °C 0 to 90% humidity (non-condensing)
eDAQ reserves the right to alter these specifications at any time.	