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## Acknowledgements

The authors would like to thank Dr. J. Kavanagh and Mr. S. Hasein for discussions, intellectual and technical input to the study of membrane fouling using electrical impedance spectroscopy.

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1. Tay, K.G. and Song, L., (2005) A more effective method for fouling characterization in a full-scale reverse osmosis process. *Desalination*, 177, 95-107.
2. Madaeni, S, Fane, A. G. and Wiley, D. (1999) Factors influencing critical flux in filtration of activated sludge. *J. Chem Tech Biotech* 74; 539-543.

3. Li, H., A.G. Fane, H.G.L. Coster and S. Vigneswaran. (1998), Direct observation of particle deposition on the membrane surface during crossflow microfiltration, *Journal of Membrane Science*. 149 pp. 83-97.
4. Zhang, Y.P. Fane, A.G. and Law, A.W.K. (2006) Critical flux and particle deposition of bidisperse suspensions during crossflow microfiltration, *Journal of Membrane Science*. 282; 189-197
5. Kwon, D.Y., Vigneswaran, S. Fane A.G. and Ben Aim R., (2000), Experimental determination of critical flux in cross-flow microfiltration, *Separation and Purification. Technology*, 19 pp. 169-181
6. Coster, H.G.L., Kim, K.J., Dahlan, K., Smith J.R., and Fell, C.J.D., (1992) Characterisation of ultrafiltration membranes by impedance spectroscopy. I. Determination of the separate electrical parameters and porosity of the skin and sublayers *Journal of Membrane Science*, 66, 19-26
7. Park, J.-S., Chilcott, T.C., Coster, H.G.L., and Moon, S.-H, (2005). Characterization of BSA-fouling of ion-exchange membrane systems using a subtraction technique for lumped data, *Journal of Membrane Science*, 246, 137-144
8. Chilcott, T.C., Chan, M., Gaedt, L., Nantawisarakul, T., Fane, A.G., and Coster H.G.L., (2002). Electrical impedance spectroscopy characterization of conducting membranes I. Theory *Journal of Membrane Science*, 195, 153-167
9. Gaedt, L., Chilcott, T.C., Chan, M., Nantawisarakul, T., Fane, A.G., and Coster. H.G.L., (2002) Electrical impedance spectroscopy characterization of conducting membranes II. *Experimental Journal of Membrane Science*, 195, 169-180
10. Kavanagh J.M., Hussein, S., Chilcott, T. C. and Coster, H. G. L. (2009). Monitoring Fouling of Reverse Osmosis Membranes using Electrical Impedance Spectroscopy *Desalination* 236; 187-193
11. Coster, H.G.L. and Chilcott, T.C. (2007a) A system for complex impedance measurement INPHAZE PCT 2006903255

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Copyright © eDAQ 2009